

## MLF -40 Freezer

The MLF -40 °C Freezer provides the reliability that you need for your valuable items.

The static refrigerated shelves help to keep an even temperature throughout. Equipped with door lock and door ajar alarm, the MLF gives you the security that you need.

### Features:

- Door Lock
- Power Failure Battery Back-up for controller
- Door Ajar Alarms
- Voltage Free Connection
- Digital Display with Built in Alarms
- Hi/Lo Read out



-40°C



115 L



600w x 600d x 830h



Model	MLF -40
Capacity	115 L
External Dimensions (including castors)	600w x 600d x 830h
Internal Dimensions	480w x 410d x 660h
Weight	100 Kg
Voltage	240V 50Hz
Watts	320w
Refrigerant	R404a
Warranty	12 Months *
Shelves size	500 x 500
Number of Shelves	3 Fixed
Digital Display	✓
Castors	X
Interior/Exterior Cabinet Finish	White Enamel
Insulation	60 mm
High/Low Temperature Alarm	✓
Static Refrigeration System	✓
Door Lock	✓
Solid Door	✓
Independent data logger	✓
Door ajar alarm	✓

\*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.