


Blanket Warmers

Features:

- Fan forced system
- TGA Approved
- Overheating alarm
- Double glazed doors
- Stainless int./ext. with chrome plated shelves
- Micro-processor control
- Temperature lock out
- Digital Display with Built In Alarms
- Hi / Lo Read out



60 °C



80 L

257 L

365 L



500w x 575d x 600h

850w x 575d x 900h

750w x 700d x 1800h



| Model | BW5 | BW10 | BW21 |
|----------------------------------|--------------------|--------------------|---------------------|
| Capacity | 80 L | 257 L | 365 L |
| External Dimensions | 500w x 575d x 600h | 850w x 575d x 900h | 750w x 700d x 1800h |
| Internal Dimensions | 390w x 480d x 430h | 740w x 485d x 590h | 630w x 610d x 1590h |
| Weight | 75. Kg | 125 Kg | 175 Kg |
| Voltage | 240 V | 240 V | 240 V |
| Warranty | 12 Months* | 12 Months* | 12 Months* |
| Shelves size | 390 x 480 | 740 x 485 | 630 x 610 |
| Number of Shelves | 2 | 3 | 5 |
| Digital Display | ✓ | ✓ | ✓ |
| Castors | X | ✓ | ✓ |
| Interior/Exterior Cabinet Finish | Stainless Steel | Stainless Steel | Stainless Steel |
| High/Low Temperature Alarm | ✓ | ✓ | ✓ |
| Glass Door | ✓ | ✓ | ✓ |

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.