


Drying Ovens

Features:

- Forced Air Convection System
- Polished Stainless steel Chamber
- Latest PID Controller
- Uniform Distribution of Air Temperature
- Double layer glass door, Larger viewing window.
- Digital Display with Built In Alarms
- Hi / Lo Read out



10°C to 250°C



80 L

136 L



740w x 618d x 630h

840w x 670d x 730h

Model	DO-9070	DO-9140
Capacity	80 L	136 L
External Dimensions	740w x 618d x 630h	840w x 670d x 730h
Internal Dimensions	450w x 400d x 450h	550w x 450d x 550h
Voltage	220 V 50 Hz	220 V 50 Hz
Power Consumption	1150 W	2050 W
Timing Range	1~9999 min	1~9999 min
Number of Shelves	2	2
Digital Display	✓	✓
Temperature Range	RT +10-200°C / RT +10-250°C	
Temperature Stability	±1.0°C	
Ambient Temperature	+5~40°C	
Warranty	12 Months*	
Optional Equipment		
Intelligent Temp Controller	✓	✓
Independent Limiting Alarm	✓	✓
R485 Connector	✓	✓
Printer	✓	✓

*Refer to warranty

Nuline Refrigeration reserves the right to alter the specifications without prior notice.