


Fluid Warmers

Features:

- Fan forced system
- TGA Approved
- Overheating alarm
- Double glazed doors
- Stainless int./ext. with chrome plated shelves
- Micro-processor control
- Temperature lock out
- Digital Display with Built In Alarms
- Hi / Lo Read out



35 °C



80 L

257 L

365 L



500w x 575d x 600h

850w x 575d x 900h

750w x 700d x 1800h



Model	FW5	FW10	FW21
Capacity	80 L	257 L	365 L
External Dimensions	500w x 575d x 600h	850w x 575d x 900h	750w x 700d x 1800h
Internal Dimensions	390w x 480d x 430h	740w x 485d x 590h	630w x 610d x 1590h
Weight	75. Kg	125 Kg	175 Kg
Voltage	240 V	240 V	240 V
Warranty	12 Months*	12 Months*	12 Months*
Shelves size	390 x 480	740 x 485	630 x 610
Number of Shelves	2	3	5
Digital Display	✓	✓	✓
Castors	X	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel	Stainless Steel	Stainless Steel
High/Low Temperature Alarm	✓	✓	✓
Glass Door	✓	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.