


MD -86 Freezer

Features:

- Microprocessor controlled temperature
- Equipped with inner doors
- Door lock
- Sensor error
- Digital Display with built in Alarms
- Hi/Lo Read out



-86 °C



588 L

688 L

838 L

936 L



872w x 1192d x 1995h

1354w x 1072d x 1907h

1160w x 1172d x 1907h

1240w x 1000d x 2000h

Optional equipment

Temperature

Data Logger



Keypad

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

Model	MD600	MD700	MD800	MD900
Capacity	588 L	688 L	838 L	936 L
External Dimensions	872w x 1192d x 1995h	1354w x 1072d x 1907h	1160w x 1172d x 1907h	1240w x 1000d x 2000h
Internal Dimensions	606w x 738d x 1310h	990w x 607d x 1140h	894w x 718d x 1310h	1000w x 700d x 1335h
Inner Doors	2	2	2	4
Weight	275 Kg	290 Kg	334 Kg	330 Kg
Watts	800 W	820 W	840 W	1600 W
Power Supply	220V/50~60 Hz;; 115V/60 Hz	220V/50~60 Hz;; 115V/60 Hz	220V/50~60 Hz;; 115V/60 Hz	220V/50~60 Hz;; 115V/60 Hz
Refrigerant	R190 & R170	HFC	R190 & R170	HFC
Ambient Temperature	32 °C	32 °C	32 °C	32 °C
Warranty	12 Months*	12 Months*	12 Months*	12 Months*
Number of Shelves	3 Adjustable	3 Adjustable	3 Adjustable	4 Adjustable
Controller	Microprocessor	Microprocessor	Microprocessor	Microprocessor
Digital Display	✓	✓	✓	✓
Castors	✓	✓	✓	✓
Hi/Lo Alarms	✓	✓	✓	✓
Power Failure Alarm	✓	✓	✓	✓
Voltage Abnormal	✓	✓	✓	✓
Sensor Error	✓	✓	✓	✓
Automatic Clean - Filter Alarm	✓	✓	✓	✓
Controller Error	✓	✓	✓	✓
High Ambient Temperature Alarm	✓	✓	✓	✓
Door Ajar Alarm	✓	✓	✓	✓
Port Entry	1 x 25 mm	2 x 25 mm	1 x 25 mm	2 x 18 mm
USB Port	✓	✓	✓	✓
Remote Alarm Port	✓	✓	✓	✓
Noise Level	55dB	55dB	55dB	55dB