

MF BTS Static Laboratory Freezers

Features:

- Top Mounted System
- TGA Approved
- Cabinet Finish In Stainless Steel Interior/Exterior
- Automatic Defrost And Water Disposal
- 5 Adjustable Shelves Per Door
- Digital Display with Built In Alarms
- Hi / Lo Read Out



TGA
Health Safety
Regulation



-18°C to -22 °C




537 L
1173 L



680w x 845d x 2000h
1340w x 845d x 2000h

Optional equipment



Model	MF600BTS	MF1200BTS
Capacity	537 L	1173 L
External Dimensions (with castors)	680w x 845d x 2000h	1340w x 845d x 2000h
Internal Dimensions	560w x 680d x 1380h	1220w x 680d x 1380h
Weight	74 Kg	169 Kg
Voltage	240 V - 50 Hz	240 V - 50 Hz
Power Requirement	10 amp GPO	10 amp GPO
Watts	650W	1250W
Refrigerant	GAS R404A / R507	GAS R404A / R507
Warranty	12 Months*	12 Months*
Shelves size	640 x 535	640 x 535
Number of Shelves	5 Adjustable	10 Adjustable
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel	Stainless Steel
Insulation	60 mm	60 mm
High/Low Temperature Alarm	✓	✓
Static (Fan Assisted) Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Tray Slides	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.