

MLS 125 Spark Safe Refrigerator

Features:

- Static Refrigeration System
- TGA Approved
- Voltage Free Connection
- Door Lock
- Door Ajar Alarms
- Automatic Condensate water Disposal
- Battery backup for controller
- Digital Display with Built In Alarms
- Hi / Lo Read out



2°C to 8°C



125 L



595w x 595d x 820h

Optional equipment



| | |
|---|---------------------------|
| Model | MLS125 |
| Capacity | 125 L |
| External Dimensions (including castors) | 595w x 595d x 820h |
| Internal Dimensions | 510w x 410d 680h |
| Weight | 42 Kg |
| Voltage | 240V 50Hz |
| Watts / AMPS | 135w / 1.05A |
| Refrigerant | R134a |
| Warranty | 12 Months* |
| Shelve Size | (3) 495x405 & (1) 495x260 |
| Number of Shelves | 3 Adjustable |
| Digital Display | ✓ |
| Interior/Exterior Cabinet Finish | White Enamel |
| Insulation | 44mm |
| High/Low Temperature Alarm | ✓ |
| Static Refrigeration System | ✓ |
| Door Lock | ✓ |
| Solid Door / Glass Door | ✓ |
| Independent data logger | ✓ |
| Optional Equipment | |
| Door ajar alarm | ✓ |
| Port entry | ✓ |
| Remote alarm | ✓ |
| Temp alert | ✓ |
| Chart recorder | ✓ |
| Battery backup for controller | ✓ |
| Castors | ✓ |

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.