

## MLi 125 Refrigerator Incubator

### Features:

- Fan Forced System
- TGA Approved
- Voltage Free Connection
- Door Lock
- Door Ajar Alarms
- Interior LED Lighting
- Automatic Condensate water Disposal
- Fan Door Switch
- Data Logger Included
- Battery Backup for Controller
- Digital Display with Built In Alarms
- Hi / Lo Read out



10°C to 40°C




125 L



595w x 595d x 820h

## Optional equipment



Temperature

Data Logger

Keypad

<b>Model</b>	<b>MLi125</b>
Capacity	125 L
External Dimensions (including castors)	595w x 595d x 820h
Internal Dimensions	510w x 410d 680h
Weight	42 Kg
Voltage	240V 50Hz
Watts	135w / 1.05A
Refrigerant	R134a
Warranty	12 Months*
Shelve Size	(3) 495x405 & (1) 495x260
Number of Shelves	3 Adjustable
Digital Display	✓
Interior/Exterior Cabinet Finish	White Enamel
Insulation	44mm
High/Low Temperature Alarm	✓
Fan Forced Refrigeration System	✓
Door Lock	✓
Solid Door / Glass Door	✓
Independent data logger	✓
<b>Optional Equipment</b>	
Door ajar alarm	✓
Port entry	✓
Remote alarm	✓
Temp alert	✓
Chart recorder	✓
Battery backup for controller	✓
Castors	✓

\*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.