


NHF Spark Safe Freezer

Features:

- Top Mounted System
- TGA Approved
- Heated Cabinet facia to prevent condensation
- Auto Defrost
- Castors / Legs
- Stainless Steel Interior
- Digital Display with Built-in Alarms
- Adjustable Shelves
- Hi / Lo Read Out



-18 °C to -22 °C



356 L


505 L



650w x 700d x 1855h

800w x 750d x 1900h

Optional equipment



Temperature

Data Logger

Keypad

Model	NHF400	NHF600
Capacity	356 L	505 L
External Dimensions (without castors)	650w x 700d x 1855h	800w x 750d x 1900h
Internal Dimensions	470w x 535d	590w x 535d
Weight	95 Kg	115 Kg
Voltage	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO
Watts	780W	960W
Refrigerant	R404A	R404A
Warranty	12 Months*	12 Months*
Shelves size	450 x 515	570 x 515
Number of Shelves	5	5
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel / White Colourbond	Stainless Steel / White Colourbond
Insulation	75 mm	75 mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.