

NHFTS Top Mounted Spark Safe Freezer

Features:

- Top Mounted System
- TGA Approved
- Heated Cabinet facia to prevent condensation
- Auto Defrost
- Castors / Legs
- Stainless Steel Interior
- Digital Display with Built-in Alarms
- Hi / Lo Read Out



-18 °C to -22°C



350 L

570L



650w x 750d x 1915h

770w x 750d x 1965h

Optional equipment



Temperature

Data Logger

Keypad

Model	NHFTS400	NHFTS600
Capacity	350 L	570 L
External Dimensions (including castors)	650w x 750d x 1915h	770w x 750d x 1965h
Internal Dimensions	470w x 535d x 1405h	590w x 535d x 1455h
Weight	110 Kg	130 Kg
Voltage	240V 50Hz	240V 50Hz
Watts	720w	960w
Power Requirement	10amp GPO	10amp GPO
Power Supply	3 Pin plug	3 Pin plug
Refrigerant	R404a	R404a
Warranty	12 Months *Refer to warranty sheet	12 Months *Refer to warranty sheet
Shelves size	450 x 515	570 x 515
Number of Shelves	5 Adjustable	5 Adjustable
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	Int Stainless Steel / Ext. White Colourbond	Int Stainless Steel / Ext. White Colourbond
Insulation	75mm	75mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.