


NHRS Spark Safe Refrigerator

Features:

- Bottom Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Glass or Solid Door
- Auto Water Disposal
- Castors
- Door Lock
- Digital Display with built in Alarms
- Hi / Lo Read out



2°C to 8°C



350 L

590 L



600w x 620d x 1850h

770w x 750d x 1900h

Optional equipment



Temperature

Data Logger

Keypad

Model	NHRS400	NHRS600
Capacity	350 L	590 L
External Dimensions (without castors)	600w x 620d x 1850h	770w x 750d x 1900h
Internal Dimensions	500w x 470d x 1750h	670w x 600d x 1740h
Weight	75 Kg	95 Kg
Voltage	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO
Watts	480W	575W
Refrigerant	R134a	R134a
Warranty	12 Months*	12 Months*
Shelves size	480 x 450	650 x 580
Number of Shelves	5	5
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond
Insulation	50 mm	50 mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door / Glass Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.