

NHRTS Top Mounted Spark Safe Refrigerator

Features:

- Digital Display with built in Alarms
- Adjustable Shelves
- Door Lock
- Hi/Lo Read out
- E-Glass & Heated cabinet facia to prevent condensation



2 °C to 8 °C



350 L


590 L



600w x 620d x 1915h

770w x 750d x 1965h

Optional equipment



Temperature

Data Logger

Keypad

Model	NHRTS400	NHRTS600
Capacity	350 L	590 L
External Dimensions (including castors)	600w x 620d x 1915h	770w x 750d x 1965h
Internal Dimensions	500w x 470d x 1490h	670w x 600d x 1530h
Weight	95 Kg	115 Kg
Voltage	240V 50Hz	240V 50Hz
Watts	386w	480w
Power Requirement	10amp GPO	10amp GPO
Refrigerant	R134a	R134a
Warranty	12 Months*	12 Months*
Shelves size	480 x 450	650 x 580
Number of Shelves	5 Adjustable	5 Adjustable
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond
Insulation	50mm	50mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door / Glass Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.