

NLD Laboratory Display Freezer

Features:

- Bottom Mounted System
- TGA Approved
- Upright Display Laboratory Freezer
- Automatic Defrost and Water Disposal
- Interior Lighting
- Castors
- Door Locks
- Digital Display with Built in Alarms
- Hi / Lo Read out



-18°C to -22 °C



412 L


920 L



680w x 720d x 1990h

1370w x 720d x 1990h

Optional equipment



Model	NLDF420	NLDF930
Capacity	412 L	920 L
External Dimensions	680w x 720d x 1990h	1370w x 720d x 1990h
Internal Dimensions	560w x 445d x 1330h	1250w x 445 d x 1330h
Weight	95 Kg	150 Kg
Voltage	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO
Watts	1000W	1100W
Refrigerant	R404a	R404a
Warranty	12 Months*	12 Months*
Shelves size	550 x 490	615 x 445
Number of Shelves	5	10
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	White Enamel	White Enamel
Insulation	60 mm	60 mm
High/Low Temperature Alarm	✓	✓
Fan Forced Refrigeration System	✓	✓
Door Lock	✓	✓
Glass Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.