


# NLMB Laboratory and Medical

## Features:

- Bottom Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Mullion Lighting
- 4 P.V.C Adjustable full width Shelves
- Double Glazed Doors
- Automatic Condensate water Disposal
- Digital Display with Built In Alarms
- Hi / Lo Read out



2°C to 8°C



488 L  
820 L  
1171 L



870w x 700d x 1980h  
1200w x 700d x 1980h  
1880w x 700d x 1980h

## Optional equipment



Model	NLMB/1	NLMB/2	NLMB/3
Capacity	488 L	820 L	1171 L
External Dimensions	750w x 700d x 1885h	1200w x 700d x 1885h	1660w x 700d x 1885h
Internal Dimensions	770w x 600d x 1450h	1100w x 600d x 1450h	1780w x 600d x 1450h
Weight	95 Kg	125 Kg	175 Kg
Voltage	240 V	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO	10 amp GPO
Watts	775W	984W	1012W
Refrigerant	R134a	R134a	R134a
Warranty	12 Months*	12 Months*	12 Months*
Shelves size	750 x 530	1080 x 530	585 x 530
Number of Shelves	4	4	12
Digital Display	✓	✓	✓
Castors	✓	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond	White Colourbond
Insulation	50 mm	50 mm	50 mm
High/Low Temperature Alarm	✓	✓	✓
Fan Forced Refrigeration System	✓	✓	✓
Door Lock	✓	✓	✓
Solid Door / Glass Door	✓	✓	✓
<b>Optional Equipment</b>			
Door ajar alarm	✓	✓	✓
Port entry	✓	✓	✓
Remote alarm	✓	✓	✓
Temp alert	✓	✓	✓
Independent data logger	✓	✓	✓
Chart recorder	✓	✓	✓
Battery backup for controller	✓	✓	✓
Additional shelves	✓	✓	✓

\*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.