

## NLMBS Spark Safe Refrigerator

### Features:

- Bottom Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Glass or Solid Door
- Auto Water Disposal
- Castors
- Door Lock
- Digital Display with built in Alarms
- Hi / Lo Read out



2°C to 8 °C



820 L

1171 L




1200w x 700d x 1980h

1880w x 700d x 1980h



## Optional equipment



Model	NLMBS/2	NLMBS/3
Capacity	820 L	1171 L
External Dimensions (without castors)	1200w x 700d x 1980h	1880w x 700d x 1980h
Internal Dimensions	1100w x 600d x 1450h	1780w x 600d x 1450h
Weight	125 Kg	175 Kg
Voltage	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO
Watts	984W	1012W
Refrigerant	R134a	R134a
Warranty	12 Months*	12 Months*
Shelves size	1050 x 530	585 x 530
Number of Shelves	4	12
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond
Insulation	50 mm	50 mm
High/Low Temperature Alarm	✓	✓
Static Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door / Glass Door	✓	✓
<b>Optional Equipment</b>		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Additional shelves	✓	✓

\*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.