


NLMS Spark Safe Laboratory Refrigerator

Features:

- Top Mounted System
- TGA Approved
- E-Glass & Heated cabinet facia to prevent condensation up to 80% R.H.
- Glass or Solid Door
- Auto Water Disposal
- Castors
- Door Lock
- Digital Display with built in Alarms
- Hi / Lo Read out



2°C to 8 °C



700 L

1000 L

1614 L

2170 L



870w x 700d x 2100h

1200w x 700d x 2100h

1880w x 700d x 2100h

2400w x 700d x 1980h

Optional equipment



Temperature

Keypad

Model	NLMS 700/1	NLMS 1000/2	NLMS 1614 /3	NLMS 2170 /4
Capacity	700 L	1000 L	1614 L	2170 L
External Dimensions (without castors)	870w x 700d x 2100h	1200w x 700d x 2100h	1880w x 700d x 2100h	2400w x 700d x 1980h
Internal Dimensions	770w x 550d x 1555h	1100w x 550d x 1555h	1780w x 550d x 1520h	2300w x 550d x 1520h
Weight	95 Kg	125 Kg	175 Kg	195 Kg
Voltage	240V 50Hz	240V 50Hz	240V 50Hz	240V 50Hz
Power Requirement	10amp GPO	10amp GPO	10amp GPO	10amp GPO
Watts	850W	1000W	1128W	1970W
Refrigerant	R134a	R134a	R134a	R134a
Warranty	12 Months*	12 Months*	12 Months*	12 Months*
Shelves size	750 x 530	1080 x 530	565 x 530	575 x 530
Number of Shelves	5 Adjustable	5 Adjustable	15 Adjustable	20 Adjustable
Digital Display	✓	✓	✓	✓
Castors	✓	✓	✓	✓
Interior/Exterior Cabinet Finish	White Colourbond	White Colourbond	White Colourbond	White Colourbond
Insulation	50mm	50mm	50mm	50mm
High/Low Temperature Alarm	✓	✓	✓	✓
Static Refrigeration System	✓	✓	✓	✓
Door Lock	✓	✓	✓	✓
Solid Door / Glass Door	✓	✓	✓	✓
Optional Equipment				
Door ajar alarm	✓	✓	✓	✓
Port entry	✓	✓	✓	✓
Remote alarm	✓	✓	✓	✓
Temp alert	✓	✓	✓	✓
Independent data logger	✓	✓	✓	✓
Chart recorder	✓	✓	✓	✓
Battery backup for controller	✓	✓	✓	✓
Additional shelves	✓	✓	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.