

Series 800 Incubators

Features:

- Fan Forced System
- Microprocessor controller (with timing function)
- With inner glass door for easy observation
- Polished stainless steel chamber
- Over Temp Protection
- Digital Display with Built In Alarms
- Hi / Lo Read out



5°C to 65 °C



35 L

50 L

80 L



620w x 490d x 490h

690w x 500d x 500h

780w x 530d x 560h

Specifications

Model	i-132B	i-152B	i-182B
Capacity	35 L	50 L	80 L
External Dimensions	620w x 490d x 490h	690w x 500d x 500h	780w x 530d x 560h
Internal Dimensions	340w x 320d x 320h	415w x 360d x 355h	500w x 400d x 650h
Watts	200W	250W	350W
Voltage	220V - 50Hz	220V - 50Hz	220V - 50Hz
Timing Range	1 ~ 9999min	1 ~ 9999min	1 ~ 9999min
Number of Shelves	2	2	2
Digital Display	✓	✓	✓
Temperature Range		RT +5~65°C	
Display Resolution		0.1°C / ±0.5°C	
Ambient Temperature		+5~35°C	
Warranty		12 Months*	
Optional Equipment			
Intelligent Programmable Temperature Controller	✓	✓	✓
Independent Limiting Alarm	✓	✓	✓
R485 Connector	✓	✓	✓
Printer	✓	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.