


## Series 900 Incubators

### Features:

- Fan Forced System
- Auto controller of fan speed
- Shelf Check Function- easy to identify problems
- Over Temperature and temperature different alarms
- Programmable controller.
- Digital Display with Built In Alarms
- Hi / Lo Read out



5°C to 80 °C



45 L

80 L



350w x 370d x 350h

400w x 405d x 500h

Model	i9042	i9082
Capacity	45 L	80 L
Internal Dimensions	350w x 370d x 350h	400w x 405d x 500h
Voltage	220V - 50Hz	220V - 50Hz
Power Consumption	250W	300W
Timing Range	1 ~ 9999min	1 ~ 9999min
Number of Shelves	2	2
Digital Display	✓	✓
Temperature Range	RT +5~80°C	
Temperature Stability	±0.3°C	
Ambient Temperature	+5~35°C	
Display Resolution	0.1°C/±0.2°C	
Warranty	12 Months*	
Optional Equipment		
Intelligent Programmable Temperature Controller	✓	✓
Independent Limiting Alarm	✓	✓
R485 Connector	✓	✓
Printer	✓	✓

\*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.