

DW -40 °C Freezer

Features:

- Bottom Mounted System
- Cabinet Finish in Galvanised steel/powder coated
- 3-5 Adjustable Shelves Per Door depending on the capacity of the freezer
- Digital Display with Built In Alarms
- Hi / Lo Read Out



-10°C to -40 °C



280 L

360 L

480 L

590 L



620w x 650d x 1920h

620w x 780d x 1920h

780w x 780d x 1920h

920w x 750d x 1920h

Optional equipment

Temperature

Data Logger



Keypad

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

Model	DW300	DW400	DW500	DW600
Capacity	280 L	360 L	480 L	590 L
External Dimensions (including castors)	620w x 650d x 1920h	620w x 780d x 1920h	780w x 780d x 1920h	920w x 750d x 1920h
Internal Dimensions	440w x 450d x 1310h	440w x 580d x 1310h	600w x 580d x 1310h	780w x 600d x 1310h
Voltage	240V 50Hz	240V 50Hz	240V 50Hz	240V 50Hz
Watts	960w	960w	1000w	1248w
Refrigerant	R404a	R404a	R404a	R404a
Warranty	12 Months*	12 Months*	12 Months*	12 Months*
Shelves size	430 x 410	430 x 550	580 x 550	760 x 570
Number of Shelves	3	4	5	5
Inner Door	2	2	2	2
Controller	Microprocessor	Microprocessor	Microprocessor	Microprocessor
Display	LCD Touchscreen	LCD Touchscreen	LCD Touchscreen	LCD Touchscreen
Castors	✓	✓	✓	✓
Interior/Exterior Cabinet Finish	Galvanised Steel / Powder	Galvanised Steel / Powder	Galvanised Steel / Powder	Galvanised Steel / Powder
Insulation	90mm	90mm	90mm	90mm
High/Low Temperature Alarm	✓	✓	✓	✓
Static Refrigeration System	✓	✓	✓	✓
Sensor Error	✓	✓	✓	✓
Hot Condenser	✓	✓	✓	✓
Low Battery Alarm	✓	✓	✓	✓
High Ambient Temperature Alarm	✓	✓	✓	✓
Door Lock	✓	✓	✓	✓
Solid Door	✓	✓	✓	✓
Door ajar alarm	✓	✓	✓	✓
Port entry	✓	✓	✓	✓
Voltage Free Contacts	✓	✓	✓	✓
Temp alert	✓	✓	✓	✓
Built in data logger via Touchscreen	✓	✓	✓	✓