


HF Series Spark Safe Freezer

Features:

- Audible Alarm
- Static Refrigerated Shelves
- Door Lock
- Digital Display with built in Alarms
- Bottom Mounted System
- TGA Approved
- HI/Lo Read out



-18 °C to -22 °C



129 L

361 L

620 L



600w x 585d x 855h

600w x 585d x 1855h

777w x 695d x 1895h



Optional equipment



Temperature

Data Logger

Keypad

Model	HF200	HF400	HF600
Capacity	129 L	361 L	620 L
External Dimensions	600w x 585d x 855h	600w x 585d x 1855h	777w x 695d x 1895h
Internal Dimensions	510w x 485d x 620h	510w x 485d x 1620h	657w x 580d x 1660h
Weight	55 Kg	75 Kg	95 Kg
Voltage	240 V	240 V	240 V
Power Requirement	10 amp GPO	10 amp GPO	10 amp GPO
Refrigerant	R404a	R404a	R404a
Warranty	12 Months*	12 Months*	12 Months*
Shelves size	460 x 440	460 x 440	650 x 525
Number of Shelves	3	6	6
Digital Display	✓	✓	✓
Interior/Exterior Cabinet Finish	White enamel	White enamel	White enamel
Insulation	60 mm	60 mm	60 mm
High/Low Temperature Alarm	✓	✓	✓
Static Refrigeration System	✓	✓	✓
Door Lock	✓	✓	✓
Solid Door	✓	✓	✓
Optional Equipment			
Door ajar alarm	✓	✓	✓
Port entry	✓	✓	✓
Remote alarm	✓	✓	✓
Temp alert	✓	✓	✓
Independent data logger	✓	✓	✓
Chart recorder	✓	✓	✓
Battery backup for controller	✓	✓	✓
Castors	-	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.