


MF Series - Laboratory Freezers

Features:

- Top Mounted System
- TGA Approved
- Cabinet Finish In Stainless Steel Interior/Exterior
- Automatic Defrost And Water Disposal
- 5 Adjustable Shelves Per Door
- Digital Display with Built In Alarms
- Hi / Lo Read Out



-18°C to -22 °C



700 L

1400L



740w x 830d x 2010h

1480w x 830d x 2010h

Optional equipment



Model	MF70BT	MF140BT
Capacity	700 L	1400 L
External Dimensions (without castors)	740w x 830d x 2010h	1480w x 830d x 2010h
Internal Dimensions	585w x 650d x 1430h	1340w x 650d x 1400h
Weight	122 Kg	159 Kg
Voltage	240 V	240 V
Watts	650 W	700 W
Power Requirement	10 amp GPO	10 amp GPO
Refrigerant	GAS R404A / R507	GAS R404A / R507
Warranty	12 Months*	12 Months*
Shelves size	530 x 640	530 x 640
Number of Shelves	5 Adjustable	10 Adjustable
Digital Display	✓	✓
Castors	✓	✓
Interior/Exterior Cabinet Finish	Stainless Steel	Stainless Steel
Insulation	70 mm	70 mm
High/Low Temperature Alarm	✓	✓
Fan Forced Refrigeration System	✓	✓
Door Lock	✓	✓
Solid Door	✓	✓
Optional Equipment		
Door ajar alarm	✓	✓
Port entry	✓	✓
Remote alarm	✓	✓
Temp alert	✓	✓
Independent data logger	✓	✓
Chart recorder	✓	✓
Battery backup for controller	✓	✓
Tray Slides	✓	✓
Additional shelves	✓	✓

*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.