


# MLB Breast Milk Refrigerator

Features:

- Fan Forced System
- TGA Approved
- Voltage Free Connection
- Door Lock
- Door Ajar Alarms
- Interior LED Lighting
- Automatic Condensate water Disposal
- Fan Door Switch
- Data Logger Included
- Battery Backup for Controller
- Digital Display with Built In Alarms
- Hi / Lo Read out



2°C to 4°C




29 L  
59 L  
125 L



445w x 485d x 535h  
445w x 485d x 725h  
595w x 595d x 825h



## Optional equipment



\*Refer to warranty sheet

Nuline Refrigeration reserves the right to alter the specifications without prior notice.

Model	MLB29GP	MLB59GP	MLB125GP
Capacity	29 L	59 L	125 L
External Dimensions	445w x 485d x 535h	445w x 485d x 725h	595w x 595d x 825h
Internal Dimensions	360w x 200d x 400h	360w x 200d x 590h	510w x 410d x 680h
Weight	26 Kg	29 Kg	38 Kg
Voltage	240V 50Hz	240V 50Hz	240V 50Hz
Watts	100w / 0.8A	100w / 0.8A	135w / 1.05A
Refrigerant	134a	134a	134a
Warranty	12 Months*	12 Months*	12 Months*
Shelves size	348 x 192	(2) 348x328 & (1) 348x192	(3) 495x405 & (1) 495x260
Number of Shelves	2 Adjustable	3 Adjustable	4 Adjustable
Digital Display	✓	✓	✓
Castors	Optional	Optional	Optional
Interior/Exterior Cabinet Finish	White Enamel	White Enamel	White Enamel
Insulation	44mm	44mm	44mm
High/Low Temperature Alarm	✓	✓	✓
Fan Forced Refrigeration System	✓	✓	✓
Door Lock	✓	✓	✓
Solid Door / Glass Door	✓	✓	✓
Remote alarm	✓	✓	✓
Independent data logger	✓	✓	✓
Chart recorder	✓	✓	✓
Additional shelves	✓	✓	✓
<b>Optional Equipment</b>			
Port entry	✓	✓	✓
Additional Shelves	✓	✓	✓
Castors	✓	✓	✓